

CLEARFIL MAJESTY™ ES-2 Universal



Intelligent Simplification

How many different shades of composite do you need to create appealing restorations in virtually every clinical situation? If you opt for CLEARFIL MAJESTY™ ES-2 Universal shade concept, a few shades will do the trick.

1

1 Shade for posterior restorations

A single universal shade is the go-to solution for all posterior restorations, independent of the colour of the underlying and adjacent tooth structure.

2

2 Shades for anterior restorations

In the aesthetically more demanding anterior region, there are two shade options to choose from.

X

No opaquer needed


A blocker or opaquer is not needed – neither in the posterior nor in the anterior region for your regular cases.

Features

- ✓ Excellent aesthetics and ease of use with a reduced number of shades and without opaquer or blocker
- ✓ Easy sculptability, non-stickiness and long working time
- ✓ Natural and long-term gloss after easy polishing
- ✓ Excellent and superior mechanical properties: low wear, low shrinkage and gentle to the opposing teeth

Benefits

- ✓ Simplified, time-saving procedure with fewer, potentially error-prone steps, for predictable results and a minimized inventory
- ✓ Convenient handling and enough time for giving the artistic touch to your restorations
- ✓ Appealing results that last through time
- ✓ Durable restorations



Excellent mechanical properties

Flexural strength	118 MPa
Filler load	78 wt%
Compressive strength	347 MPa
Volumetric shrinkage	1.9 %
Curing depth	2.0 mm
Working time under ambient light	270 sec

Source: Kuraray Noritake Dental Inc.

Shade comparison table

Universal shades are designed to cover all cavity classes with a few shades based on a following chart:

Shade category	Cavity class	Shade				
		A1 B1, B2	A2 D2, C1	A3 C2, D3	A3.5 B3, B4	A4 C3, C4
Universal	I	U Universal				
	II					
	III	UL Universal Light		UD Universal Dark		
	IV					
	V					

For clinical cases in the anterior region that require high lightness/value, we offer a dedicated white shade (UW).

Product assortment

ARTICLE NR	CLEARFIL MAJESTY™ ES-2	CONTENTS
#3181-EU	Universal Syringe (U)	1 syringe U: 2.0mL (3.6g)
#3183-EU	Universal Syringe Light (UL)	1 syringe UW: 2.0mL (3.6g)
#3184-EU	Universal Syringe Dark (UD)	1 syringe UL: 2.0mL (3.6g)
#3182-EU	Universal Syringe White (UW)	1 syringe UD: 2.0mL (3.6g)
#3191-EU	Universal PLT (U)	20 PLT U: 0.25g
#3193-EU	Universal PLT Light (UL)	20 PLT UW: 0.25g
#3194-EU	Universal PLT Dark (UD)	20 PLT UL: 0.25g
#3192-EU	Universal PLT White (UW)	20 PLT UD: 0.25g

SCAN TO LEARN MORE

