

PANAVIA™ SA Cement Universal



**Single-step
cementation,
even to glass
ceramics**

Using a unique silane coupling agent built into the paste - LCSi monomer - the cement delivers a strong, durable chemical bond to porcelain, lithium disilicate and composite resin. The original MDP monomer present in the paste too, allows for chemical reactivity with zirconia, dentin and enamel.

**A breakthrough from the innovators of adhesion
and silane technology**



The unique LCSi monomer creates a strong chemical bond to Porcelain, Lithium Disilicate and Composite Resin.



The Original MDP monomer¹ creates a strong chemical bond to Enamel, Dentin, Metal Alloy and Zirconia.

A forerunner in self-adhesive resin cements

Building on PANAVIA™ SA Cement Universal offers wide range of indications for adhesion and is suitable for crowns/bridges, inlays/onlays, posts, and even adhesion bridges. Thanks to its powerful single-step bonding formula, cementation is easier than ever before. The product is available in both automix and handmix options.



Indications

- ✓ Cementation of crowns, bridges, inlays and onlays
- ✓ Cementation of prosthetic restorations on implant abutments and frames
- ✓ Cementation of adhesion bridges and splints
- ✓ Cementation of posts and cores
- ✓ Amalgam bonding

Benefits

- ✓ Wide indication coverage including adhesion bridges and lithium disilicate
- ✓ Easy storage at room temperature*
- ✓ Easy application with little pressure
- ✓ Improved Endo tip for easier root canal application
- ✓ Automix and Handmix delivery
- ✓ Easy, gum friendly excess removal

*Below 25°C despite the presence of silane(LCSI).

Shades



Universal (A2)



White



Translucent

Product assortment

ARTICLE NR	PANAIA™ SA Cement Universal	CONTENTS
#4202-EU	Handmix Universal (A2)	1 syringe per shade (9.2g / 5.0 ml)
#4203-EU	Handmix White	1 syringe per shade (9.2g / 5.0 ml)
#4204-EU	Handmix Translucent	1 syringe per shade (9.2g / 5.0 ml)
#4206-EU	Handmix Double Pack Universal (A2)	2 syringes per shade (9.2g / 5.0 ml), Mixing spatula, Mixing pad
#4212-EU	Automix Universal (A2)	1 syringe per shade (8.2g / 4.6 ml), 20 Mixing tips
#4213-EU	Automix White	1 syringe per shade (8.2g / 4.6 ml), 20 Mixing tips
#4214-EU	Automix Translucent	1 syringe per shade (8.2g / 4.6 ml), 20 Mixing tips
#4210-EU	Automix Value Pack Universal (A2)	3 syringes per shade (8.2g / 4.6 ml), 40 Mixing tips, 10 Endo tips (SS)
#4211-EU	Automix Value Pack Translucent	3 syringes per shade (8.2g / 4.6 ml), 40 Mixing tips, 10 Endo tips (SS)
#4218-EU	Automix Mixing tip (Regular)	20 Mixing tips
#4221-EU	Automix Endo tip (SS)	20 Endo tips (SS)

1. Original MDP technology. The final result may be influenced by minor evaluation conditions stipulated by ISO regulations.

SCAN TO LEARN MORE

