











COLOR COMBINATION TABLE



EU Importer
Kuraray Europe GmbH
Philipp-Reis-Strasse 4,
65795 Hattersheim am Main, Germany
Phone +49 (0)69 305 35 835
Fax +49 (0)69 305 98 35 835
www.kuraraynoritake.eu
centralmarketing@kuraray.com



Kuraray Noritake Dental Inc.
300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
Website www.kuraraynoritake.com



09/19 MKA00584

kuraray

Noritake

Super Porcelain **EX-3**

Porcelain for metal ceramic dental restorations

Noritake Super Porcelain EX-3
For True-to-life result



Kuraray Noritake Dental Inc.

kuraraynoritake.eu

MIMIC NATURE

A TECHNOLOGICAL MASTERPIECE THAT REFLECTS THE NATURALNESS OF DENTITION THROUGH YOUR HANDS.


Noritake Super Porcelain EX-3 is a high-quality porcelain system that was developed as a result of the extensive knowledge and ceramic technology of Kuraray Noritake Dental.

PRECISION OF THE NATURAL DENTITION

Kuraray Noritake Dental's shade concept creates an intra-oral harmony by diminishing the greenish color and adding a pinkish hue.

FLUORESCENCE CREATES THE NATURAL LOOK

Noritake Super Porcelain EX-3 features a clean fluorescent pale blue which creates the look of a natural tooth.



Noritake Super Porcelain EX-3 / Natural Tooth under UV Light

CLINICAL DATA

Survival Analysis of PFM (with Noritake Super Porcelain EX-3) posterior crowns

Group	Total # Crowns	Mean Survival (SE) (Days)	Probability of Survival in %
PFM	1080	2694 (5)	99.3

Fusun Ozer, DMD, PhD others, "A retrospective survey on long-term survival of posterior zirconia and porcelain-fused-to-metal crowns in private practice", Quintessence Int 2014;45:31-28

Coping system: Lava, 3M ESPE; Procera, Nobel Biocare; KATANATM, Kuraray Noritake Dental Fabricated by one dental laboratory (Cusp, Boston)

CASE 1

Porcelain Fused to Metal Restoration



Tooth preparation



Postoperative view

Fabrication by Naoki Hayashi, RDT

CASE 2

Porcelain Laminate Veneer Restoration



Preoperative view



Postoperative view

Fabrication by Nondas Vlachopoulos, RDT

CASE 3

PFM Incisor Crowns using Noritake Super Porcelain EX-3



Preoperative view




1 PFM. 3 Laminates on refractory



Postoperative view

Fabrication by Daniele Rondoni, RDT

Noritake Super Porcelain EX-3 chromatic map



by Daniele Rondoni, RDT

COLOR COMBINATION TABLE

PASTE OPAQUE 6g / 2.4ml

POA1	POA2	POA3	POA3.5	POA4	POB1	POB2	POB3
POB4	POC1	POC2	POC3	POC4	POD2	POD3	POD4
POnA1	POnA2	POnA3	POnA3.5	POnA4	POnB1	POnB2	POnB3
POnB4	POnC1	POnC2	POnC3	POnC4	POnD2	POnD3	POnD4
PONW0	PONW0.5	PONP1.5	PONP2.5	POEW0	POEW	POBA	NP Bonder*

*When the different types of alloy is used, NP Bonder color will result in different color outcome.

PASTE OPAQUE MODIFIER 3g / 1.2ml

PO White	PO Gray	PO Blue	PO Yellow	PO Orange	PO Earth Brown	PO Reddish Brown	PO Pink
----------	---------	---------	-----------	-----------	----------------	------------------	---------

UNIVERSAL PASTE OPAQUE 6g

UPnA1	UPnA2	UPnA3	UPnA3.5	UPnA4	UPnB1	UPnB2	UPnB3
UPnB4	UPnC1	UPnC2	UPnC3	UPnC4	UPnD2	UPnD3	UPnD4
UPNW0	UPNW0.5	UPNP1.5	UPNP2.5	UPEW0	EPEW		

UNIVERSAL PASTE OPAQUE MODIFIER 3g

UP White	UP Gray	UP Blue	UP Yellow	UP Orange	UP Earth Brown	UP Reddish Brown	UP Pink
----------	---------	---------	-----------	-----------	----------------	------------------	---------

POWDER OPAQUE 10g / 50g / 200g

A10	A20	A30	A3.50	A40	B10	B20	B30
B40	C10	C20	C30	C40	D20	D30	D40
nA10	nA20	nA30	nA3.50	nA40	nB10	nB20	nB30
nB40	nC10	nC20	nC30	nC40	nD20	nD30	nD40
NW00	NW0.50	NP1.50	NP2.50	EW00	EW0		

POWDER OPAQUE MODIFIER 10g / 50g

OM Gray	OM Orange	OM Pink
---------	-----------	---------

MARGIN 10g / 50g

MA1	MA2	MA3	MA3.5	MA4	MB1	MB2	MB3
MB4	MC2	MC4	MD3	MD4	MNW0	MNW0.5	MNP1.5
MNP2.5	M Clear	M Orange	M Peach	MDL	MRP		

OPACIOUS BODY 10g / 50g 200g

OBA1	OBA2	OBA3	OBA3.5	OBA4	OBB1	OBB2	OBB3
OBB4	OBC1	OBC2	OBC3	OBC4	OBD2	OBD3	OBD4
OBNP1.5	OBNP2.5	OBWhite	OBOrange	OBPale Pink	OBE	OB1110	OB1120
OB2015	OB2025	OB2110	OB2120	OB2130	OB2215	OB2225	OB3014
OB3025	OB3110	OB3120	OB3130	OB3215	OB3225	OB4015	OB4025
OB4110	OB4120	OB4130	OB4215	OB4225	OB5110	OB5210	OB5130

CERVICAL 10g / 50g / 200g

CV-1	CV-2	CV-3	CV-4
------	------	------	------

NORITAKE SUPER PORCELAIN EX-3 LINE-UPS

Universal Paste Opaque: Same color as Paste Opaque and Paste Opaque Modifier.

Value Shade: Value Shades are additional shades that expand the existing Noritake Super Porcelain EX-3 selection of body and Opacious Body porcelain powder to include the VITA® 3D-Master® shade.

Speed Enamel (S1, S2, S3, S4): Specially developed for the two layer build-up technique.

n color Shade: n color shade is formulated to achive chroma intensified Body and Paste Opaque (Universal Paste Opaque, Powder Opaque) of the shade. With intensified chroma, it even can be used in a case where there is insufficient space for porcelain build-up.